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## Finger Stretching Trainer

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#### Instruction Manual

(SKA)

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# Contents

Product Introduction	02
User Guide	04
Safety Precautions	06
Safety Precaution Include	07
Troubleshooting	09
Technical Parameters	10

After service, Thank you for choosing the finger trainer.

Please read this manual carefully before using this product.

Please keep this manual properly for guidance when using the stretching trainer.

The finger extension trainer can maintain and expand the user's finger joint range of motion through finger circulation rehabilitation training, prevent joint contracture and muscle atrophy, and greatly promote the reorganization and compensation of the nervous system, helping patients recover as soon as possible.

The recovery of hand dysfunction after stroke is a long-term and arduous process. Patients' conditions, rehabilitation courses, and curative effects are very different. Under normal circumstances, the recovery of hand function is slower than other parts, and long-term rehabilitation training is needed. If you have any questions in use, please get in touch with [info@ziglite.com](mailto:info@ziglite.com).

If you have any questions in use, please contact: [info@ziglite.com](mailto:info@ziglite.com)

The manual version may upgrade due to changes in software version or technical specifications without prior notice.


# Product Introduction

## Scope Of Application

It is suitable for rehabilitation training in (families) institutions for patients with finger joint recovery after injury or limb dysfunction sequelae caused by brain injury, stroke, cerebral palsy, spinal cord injury, orthopedic surgery, multiple sclerosis, etc.

## Product Composition


### Trainer Host


Note: When connecting the training gloves, it is necessary to rotate and tighten the tracheal joint clockwise to the corresponding interface of the host.

# Product Introduction

## Training Gloves


## Power Adapter

# User Guide

## Buttons, Screen Icons In Detail

### Screen Display

display of training strength, contraction time, and opening time

### Function Selection Indicator Light

Set the training strength, retract the gear, and turn on the indicator light when opening the gear

### Switch Button

One click power on, pause

### Strength-

Reduce 1 intensity per click

### Strength+

Increase 1 intensity per click

### Set Switch

Press each time to switch between setting training strength, contraction time, and opening time


## Function Settings

	It can be set to 1-9 levels to meet the needs of different patients.
	1-9 levels, the higher the gear, the longer the contraction dwell time.
	1-9 levels, the higher the gear, the longer the opening dwell time.

# User Guide

## Air Valve Switch Control

Each glove of the training glove is equipped with a single independent air valve switch on the trachea, which can be opened or closed for targeted finger selection training.

### Open The Air Valve

The air valve and trachea sequence are consistent, and the air valve is open. The gas valve switch defaults to the relaxed state. The corresponding fingers can undergo training.

### Close The Air Valve

The air valve and trachea sequence are cross-shaped, and the air valve is closed; corresponding fingers are not trained.

# Safety Precautions

**[warning] Failure to comply with the following safety instructions may lead to severe accidents such as fire, electric shock, and casualties.**

- Do not disassemble, repair, or modify the components of the equipment by yourself to avoid abnormalities in the treatment device.
- Do not store components in or near water heaters, microwave ovens, cooking equipment, or high-pressure containers to avoid fire caused by overheating.
- It cannot be sterilized or autoclaved.
- Do not continue to use it when the alarm prompt appears.
- It is suitable for hand circulation and mirror movement rehabilitation training. Please do not use it for other purposes to avoid accidents.
- Do not store the parts in a place that is too hot or cold to avoid accidents.
- The validity period of the product is five years, and the disposal of the product at the end of its service life should comply with the requirements of local laws and regulations.
- Please keep the power connector dry and water-free.
- Please do not plug in or unplug the power supply with wet hands to avoid electric shock or injury.

## **Working Environment**

- 1) Ambient temperature: 10°C~35°C
- 2) Relative humidity: 30%~75%, excluding condensation
- 3) Atmospheric pressure range: 700hpa~1060hpa
- 4) Power supply: AC100~240V, 50/60Hz (trainer host power supply DC12V)

## **Storage And Transportation**

- During transportation, please fix the trainer in a relatively stable place to prevent collision and damage to the trainer.
- Please prevent liquids or foreign objects from entering the trainer to avoid damaging internal electronic components and causing malfunctions.
- Please do not place it in places with strong sunlight or ultraviolet radiation for a long time to avoid causing discoloration of the shell.
- Do not place other items on the product when saving it.
- Do not fold the gloves forcefully to avoid damaging the soft-actuated actuator and causing malfunctions.
- Do not store the trainer in a humid environment or magnetic solid field interference unsuitable for working with electronic devices.

# Safety Precaution Include

- In case of skin allergy and other discomfort while using the product, please stop using it and consult a doctor.
- In case of a sprain, trauma, bone injury of unknown condition, severe osteoporosis, untreated fracture, or people who need bracing after treatment, do not use this product without professional confirmation.
- For patients with diseases or wounds or being treated for them or with a previous medical history, please consult a doctor after becoming familiar with the following contents.

## Contraindication

- Patients with severe blood flow disorder such as thoracic venous tumors and hemocolumn disease.
- Patients with ulcer, wound, purulent fruit, bleeding, scald, inflammation, etc.

## Other Precautions

- 1) In case of abnormal experience or body feeling during use, stop using it and use it under the guidance of a doctor after seeing a doctor;
- 2) In case of redness, pruritus and rash during use, stop use and seek medical attention;
- 3) Do not use it together with other equipment or in the state of applying drugs;
- 4) Do not use it alone for unconscious people or young children;
- 5) Do not use this product for any purpose other than the purpose of use.

## Before Use

- 1) Before wearing training gloves, please try to run the therapeutic device host for 1 minute to confirm that all parts function usually;
- 2) Patients with mandatory systolic spasm need to massage before using the therapeutic device;
- 3) Note: before training, please confirm that the non muscle tension of the hand to be trained is too high (above grade II)

# Safety Precaution Include

## In Use

In case of the following conditions during operation, please stop using immediately and disconnect the power supply.

- 1) The trainer host emits abnormal sound;
- 2) The surface temperature of the trainer host is too high;
- 3) Water ingress into the trainer host;
- 4) Users experience pain, swelling, or other discomfort symptoms;
- 5) If spasms occur during training, immediately stop using it.

## Usage Time

- 1) Suggest selecting appropriate training methods based on one's actual situation;
- 2) The duration of each use should not be too long. If there is soreness in the hands the next day after training, it is necessary to reduce the number and duration of daily training;
- 3) It is recommended to use it 2-4 times a day, with a minimum interval of 45 minutes between uses;
- 4) If it is necessary to replace or extend the use time during use, please consult professional medical personnel;
- 5) Long term continuous use can easily lead to high internal temperature of the machine, which may cause phenomena such as sleep or asynchronous movement. After the machine is fully cooled, it will automatically return to normal.

## Use Of Power Adapter

- 1) The power adapter must be plugged and unplugged by a patient's supervisor;
- 2) After using the trainer, please disconnect the power adapter from the 220V power supply;
- 3) Do not use a power adapter that is not equipped with this machine.

# Troubleshooting

Problem	The Reason	Measures
Flip the "power switch" down to turn on the machine	Poor power connection	Check the power cord and power adapter, the connection between the power adapter and the host
Training gloves are not working	Gloves that are not properly attached can damage the host or glove structure	Check the physical status of the training gloves and the connection between the trainer host and the training gloves
Trainer gloves cannot be removed	Training gloves trachea joint is not rotated out	Twist out the training gloves trachea joint counterclockwise

**Note:** If you still can't use the above method, please contact us.

## Care & Maintenance

### How to clean ?

- Do not wash training gloves and mirrored gloves directly with water. Dip a soft cloth in water or neutral detergent, wring it out, and wipe the gloves.
- When dust or dirt appears on the surface of the trainer host, users can use medical alcohol with a concentration of 75% to wipe it. When wiping, please dip a small amount of alcohol into a dry cloth to avoid allowing alcohol to drip or flow into the product.
- Clean with clean and dry soft cloth, and dry in a cool and ventilated place.

### Note

- The use of cleaning and disinfecting agents other than alcohol may damage the product, shorten the service life or cause safety risks.
- It is forbidden to use high-temperature drying and exposure to the sun for quick drying or disinfection.
- Do not use gasoline, thinner or similar organic solvents to clean any part of this product. Do not use radiation, steam, ethylene oxide and other methods for disinfection, otherwise the product will be damaged.

# Technical Parameters

## Product Specifications

Product Name	Finger Stretching Trainer	Power	12W
Product Model	SKA	Force Setting	1-9
Material	Aluminium Alloy Casing	Time Setting	15 Minutes
Power Supply	AC100-240V 50/60Hz	Running Noise	0.47kg
Size	103*90*60mm	Training Mode	Circuit Training / Finger Selection Training

## Glove Specifications


Product Name	Training Gloves	Size	Left / Right S,M,L,XL,XXL
Material	PVC+Soft Silicone+Elastic Cloth	Weight	0.22kg

# Technical Parameters

## Glove Size

Size	S	M	L	XL	XXL
Hand Length (mm)	145-165	160-185	180-195	195-210	210-235
Hand Width (mm)	65-80	80-95	90-105	100-115	115-125
Middle Finger Circumference (mm)	40-50	50-60	60-70	70-90	90-100

- The corresponding model can be selected according to the actual hand measurement of the user. If two or more measured values are at the critical point, it is recommended to select one size larger.
- Women generally choose s/m, men generally choose l/xl, and the larger hand is XXL.

## Customer Service

Email : [info@ziglite.com](mailto:info@ziglite.com)

Address : Room B3, 18/F Bonsun Industrial Building, 364-366 Sha Tsui Road,  
Tsuen Wan, Hong Kong

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Office Address : Room B3, 18/F Bonsun Industrial Building, 364-366 Sha Tsui Road, Tsuen  
Wan, Hong Kong

Email : [info@ziglite.com](mailto:info@ziglite.com)

Fax: +852 3544 0462

Office : +852 3107 7500

Whatsapp : +852 6756 7703

NO. SKA03-CH-EN

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Ziglite Smart Health Care Products Co., Ltd.